

**ESS CORE N' ONE™
5 GRAM CAPSULE AND HANDLE**

- CNO is both a soil collection and transportation device
- Screw-threaded Septa closure ensures airtight seal and “zero headspace”
- CNO Capsule incorporates extruding feature; no need for expensive tools
- CNO Handle is strong, economical and easy to use

**ESS LOCK N' LOAD™
HANDLE AND SYRINGE
(For Field Preservation Technique)**

- Handle locks at 5 and 10 gram settings with one easy turn
- Beveled-edge syringe is stronger than cut syringes and fits in the neck of 40mL glass vials
- Turn handle and push to dispense soil sample into vial — no need to remove the syringe
- Strong, accurate, easy and economical

5035A SAMPLE KITS

For Field Preservation Technique

Offered in Multiple Configurations with choices of:

- Sodium Bisulfate, Methanol, UPBW DI+™ and TSP vials
- HDPE moisture container
- LOCK N' LOAD® syringe (handle sold separately)
- TERRACORE® sampler

ESS PrePreserved® VIALS

- Sodium Bisulfate, UPBW DI+™ and TSP vials for low level concentration
- Methanol for high concentration testing
- Magnetic stir bars
- Tare weight on attached label
- ESS Low Bleed Septa
- Graduation in mL available